# CERAPOWER 2.0

Instructions for use

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**invicon**  
chemical solutions

# PRODUCT INSERT – CONTENTS

|   |    |
|---|----|
| 1. Unit                                 |    |
| 1.1 Description of unit                 | 4  |
| 1.2 Technical data                      | 4  |
| 2. Instruction                          |    |
| 2.1 Safety advice                       | 5  |
| 2.2 Initiation/scope of delivery        | 6  |
| 2.3 Keyboard unit                       | 7  |
| 2.4 Operation CeraPower 2.0             | 7  |
| 2.5 Abort/termination of programme      | 8  |
| 2.6 Change between °C and °F            | 8  |
| 2.7 Reactivation of temperature monitor | 8  |
| 2.8 Care and maintenance                | 9  |
| 2.9 Exchange of fuses                   | 9  |
| 2.10 Glycerine extractor                | 9  |
| 2.11 General notes                      | 10 |
| 2.12 Guarantee / Liability              | 10 |
| 2.13 Trouble shooting                   | 11 |
| 2.14 Spare part list                    | 12 |
| 3. Explosion sketch                     | 14 |
| 4. Operation diagram                    | 15 |
| 5. Pneumatic plan                       | 16 |
| 6. Declaration of Conformity            | 17 |

# PRODUCT INSERT – DETAILS

## Definition symbols

Symbols used in this document are defined as follows.

Symbol for risk of danger:

- risk to life
- risk of physical injury
- risk of minor personal injury

Symbol for risk of danger:

- risk of property damage
- risk of environmental contamination

Symbol for useful tip or information

Symbol for action required

## Intended use

Use of this equipment in ways other than those described in this User Guide may result in injury to persons or damage to property. Use this equipment only as described in this User Guide.

invicon cannot be responsible for injuries or damages resulting from unintended applications of its equipment. Unintended uses may result from taking the following actions:

- Making changes to equipment that has not been recommended in the User Guide
- Using incompatible or damaged replacement parts
- Using unapproved accessories or auxiliary equipment

# 1. UNIT

## 1.1 Description of unit

CeraPower 2.0 units are machines for pressure-polymerisation of HyCeram which polymerise at a temperature of 120 °C under 4 bar pressure within 60 minutes.

## 1.2 Technical data

|                            |   |
|----------------------------|---|
| H x W x D                  | 230 x 290 x 330 mm                        |
| Empty weight               | 10 kg                                     |
| Power consumption          | 450 W                                     |
| Min. Operating pressure    | 3 bar                                     |
| Max. Operating pressure    | 4,2 bar                                   |
| Max. Operating temperature | 120 °C                                    |
| Volume                     | 2 Liter                                   |
| Supply voltage             | 230 V/50 Hz<br>220 V/60 Hz<br>115 V/60 Hz |

(Voltage/frequency see type plate)

|                          |        |
|--------------------------|--------|
| Min. compressor pressure | 4 bar  |
| Max. compressor pressure | 10 bar |

## 2. INSTRUCTION

### 2.1 Safety advice

Read the following notes before connecting and putting the unit into service. The operating service and the function of the unit can only be guaranteed for, if the general safety and accident prevention instructions of the legislation as well as the safety instructions of the working instruction are observed.

- This unit must only be used according to the present working instructions. We cannot be held liable for damages which have been caused by improper treatment resp. faulty manners of operation.
- The unit has to be placed onto a stable and even surface, which is suitable for the operation weight of the unit.
- The indicated voltage on the type plate has to meet the voltage of the power source.
- Connect the unit only to a socket with protective conductor. Never touch the plug with wet hands.
- The pressure of the compressed air line must not exceed 10 bar.
- The heating of the unit may only be switched on, if a sufficient quantity of glycerine has been poured.
- Do not lock the operation elements of the unit.
- Plates and labels have to be maintained in a well legible condition and must not be removed.
- Authorised operators: The operating company of the machine has to hold the working instructions at the operators' disposal and has to make sure that authorised operators have read and understood the manual. Only in that case the operators may operate the unit.
- The unit has to be examined before operation in view of proper condition. If the unit is not in proper condition, it must not be used and has to be marked accordingly.
- In case of longer works at maximum operation temperature the lid of the unit may heat up to. Therefore it must not be touched – not even for a short moment.
- At all processing steps which are executed under pressure please wait with the opening of the lid until the compressed air has completely escaped from the pressure pot.
- Please do not bend over the unit while opening the lid.
- Before cleaning and maintenance of the unit or the changing of parts the plug has to be disconnected.
- Do not intrude any objects into the unit.
- The opening of the unit and repairs are only to be carried out by accordingly trained experts.
- Arbitrary rebuilding or modifications are not permitted due to safety reasons.
- Only those appliances and spare parts may be used which are released by the manufacturer. We are not responsible for damages caused by the use of outside parts.
- If the unit is not used for a longer time, please disconnect the mains.
- According to the pressure tank regulation the CeraPower 2.0 is a pressure tank of the test group I. These pressure tanks do not need to be taken a test again. But we recommend to check the pressure pot every 3 years by the manufacturer.
- The required operation and maintenance conditions mentioned in these working instructions are to be complied with imperatively. When working with the CeraPower 2.0 the general accident prevention regulations have to be observed.

## 2. INSTRUCTION

### 2.2 Initiation

Before initiation make sure that the voltage indicated on the type plate meets the voltage of the power source.

When selecting a suitable place for the machine, please make sure that the storage space is even, resistant to splash water as well as sufficiently stable for a minimum weight of 21 kg.

Connect the unit to an air compressor. For that purpose push the pressure hose that is included in delivery onto the hose clip (25) of the unit and fix it by means of a hose clamp. Assemble the fast coupling at the other end of the tube and connect it with the pressure line.

The pressure of the compressed air line must not exceed 10 bars.

Put the power cord into the female connector (15) and create a connection with the Euro/USA plug. Open the lid (2) by means of the handle (1) and fill the unit with glycerine. (Maximum glycerine level: 30 mm below the upper edge).


Before first initiation please read the working instructions thoroughly. In case of comprehension problems contact your responsible dealer or Dreve, if necessary.

### Scope of delivery

- 1 CeraPower 2.0
- 1 power cord
- 1 air connection hose incl. hose clamps and fast coupling
- 1 glycerine drain hose incl. coupling
- 1 working instructions in 5 languages

## 2. INSTRUCTION

### 2.3 Keyboard unit

#### No. Description

- a Manometer
- b Adjustment knob
- c LED – heating
- d Temperature indication display
- e Temperature switch
- f Time indication display
- g LED – polymerisation
- h Lid opening control

### 2.4 Operation CeraPower 2.0 1/120 °C

Besides the general temperature area of 0-95 °C CeraPower 2.0 also permits to increase the glycerine temperature to 120 °C. In order to avoid that hot glycerine of 120 °C escapes from the pressure pot and in order to reduce the risk of scalding when opening the lid, the CeraPower 2.0 is equipped with a safety lock which only allows the lid to open when the operating pressure in the unit has fallen and the temperature is lower than 95 °C (203 °F). This is also applicable in case of power outage.

This means that in the application process the pressure pot has to remain closed when the unit is turned off and pressure is exerted.

When the unit is turned on at the main switch (14) the pressure pot exhausts automatically and the lid (2) can be opened.

Adjustment of temperature and polymerisation time is to be executed as described from point 7 et seqq. If a temperature > 95 °C is required, 95 °C should be selected as pre-warming temperature; the final temperature is to be adjusted directly before the polymerisation process.

According to the temperature of the glycerine bath the objects have to be inserted and taken out with suitable instruments!

The time until the final temperature is reached has to be added to the general polymerisation time. This can differ between 10 and 20 minutes, depending on the glycerine filling level. As long as the value in the temperature indication display (d) blinks, the final temperature has not been reached. The present temperature of the glycerine bath can be checked at all times by pressing the temperature switch (e), (see point 7.4).

When the time has passed the unit ventilates automatically. The lid of the CeraPower 2.0 can, however, only be opened, if the temperature of the glycerine bath has fallen below 95 °C.

## 2. INSTRUCTION

If the pressure pot has got a temperature of more than 95 °C after the time has run out, it is cooled down by cold air until the temperature has fallen below 95 °C. During this time the cooling air is released by means of the ventilation valve (26). When operating at 120 °C the cool down phase can take 5-8 minutes from time expiration to the opening of the lid.

The possibility of opening the lid is indicated by the green control lamp (h). Open the lid and remove the objects.

### 2.5 Abort/termination of programme

The polymerisation process can be terminated prematurely. For this purpose push the adjustment knob (b) for 3 seconds. The CeraPower 2.0 ventilates automatically, the display for time shows 0, the pre-adjusted temperature is maintained.

For following polymerisation processes with identical temperature continue as described from point 7.2 on.

### 2.6 Change between °C and °F

At delivery the CeraPower 2.0 units are adjusted to °C. In order to change to °F push the temperature switch (e) and the main switch (14) at the same time. The temperature indication display (d) shows L (which stands for °C). The time indication display (f) shows the maximum values 120 °C or 248 °F. Push the adjustment knob (b) one time in order to confirm the selected temperature denomination.

The time and temperature indication displays (d + f) indicate 0. Then go on as described from point 2.5 on.

### 2.7 Reactivation of temperature monitor

In order to avoid an overheating and consequently the destruction of important CeraPower 2.0 components, all units have been equipped with a temperature control (11). This part is situated from outside at the bottom of the polymerisation pot and interrupts the circuit of the heating in case of too much temperature development. In order to reactivate the temperature monitor disconnect the unit from the mains and let it cool down. Drain the glycerine completely through the drain tap (23). Lay the CeraPower 2.0 on its side and use a soft tissue as underlay, if necessary.

Unscrew the perforated metal plate at the bottom of the unit and push the red reset knob of the temperature monitor until it latches completely. Insert perforated metal plate again and tighten it. Put the CeraPower 2.0 into normal position, close drain tap (23), fill the unit with glycerine and connect to the mains by means of the power cord. Then go on as described from point 2.5 on.

## 2. INSTRUCTION

### 2.8 Care and maintenance

Disconnect the unit from mains supply before care and maintenance! The glycerine inside the pressure pot has to be changed regularly.

For the cleaning of the pressure pot open the drain tap (23) at the rear of the unit and let the glycerine flow out. Remove dirt thoroughly.


Do not use scouring agents, steel-wool or sponge scourers.

Clean the interior pot with fresh water, a neutral washing-up liquid and a soft scarf or sponge. Afterwards let the water run out. Close the drain tap (23) and fill with fresh glycerine. The O-ring should be greased regularly with silicone grease.

### 2.9 Exchange of fuses

► Step 1: Disconnect the unit from mains supply. The fuse carrier is situated at the rear of the unit, integrated in the mains supply. In order to control or exchange the fuses push the push catch of the carrier, keep it pressed and remove the carrier from the unit.


► Step 2: Pull the defective fuse out of the clamp carrier applying slight traction and replace it by a new, equivalent one.

Only use the fuse strength indicated on the rear.

► Step 3: After exchanging the fuse insert the fuse carrier into the housing applying slight pressure until the fixation latches noticeably.

### 2.10 Glycerine extractor

- Deflate the glycerine extractor solely under cold condition before operating
- Do not deflate the glycerine extractor during operating under pressure
- Glycerine extractor may heat up during operation

## 2. INSTRUCTION

### 2.11 General notes

Repairs on the CeraPower 2.0 units may only be executed by experts. Disconnect the unit from mains supply at any rate. Dreve company is only responsible for the safety, reliability and performance of the unit, if:

- the unit is exclusively used for the corresponding field of application.
- operation is only carried out by accordingly trained personnel.
- the unit is operated in accordance with these working instructions.
- extensions, new adjustments, changes or repairs are executed by persons authorised by the company.

As we are currently developing our products, we reserve the right for technical changes.

### 2.12 Guarantee

The guarantee time amounts to two years after the date of delivery. Guarantee is only accepted in case of material or fabrication defects. The guarantee obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

### Liability

The right of the purchaser to bring actions arising out of defects shall in all cases be in lapse one year after the date of delivery. The warranty obligation does not refer to natural wear and tear or to damage occurring after the transfer of risk due to improper or unsuitable handling, excessive stressing, unsuitable operating material, and chemical, electrochemical or electrical influences of a nature not provided for according to the contract. On improper alterations, repairs or maintenance or breach of seals carried out by the purchaser or third parties our warranty obligations and liability for consequences caused thereby are rendered void.

## 2. INSTRUCTION

### 2.13 Trouble shooting

#### Problem

- Unit does not show any function

- Unit does not heat

- Unit does not create pressure

- Unit loses compressed air under the lid

- Unit blows air out of the security valve

- Lid cannot be opened

#### Cause

- Mains supply is not correct
- Main switch is off
- Fuses in fuse carrier are defective

- Temperature regulator shows 0 °C
- Heating is defective

- Temperature monitor has switched off the heating

- Unit does not have compressed air
- Timer is not turned on
- Magnetic valve is defective

- Lid O-ring is polluted
- Lid O-ring is defective

- Input pressure is too high

- Temperature in the pressure tank is too high

#### Aid

- Control the mains supply
- Turn on the main switch
- Exchange fuses; in case of unsuccessful exchange look for other causes

- Adjust the temperature regulator
- Heating has to be exchanged by an expert

- Push reset button at temperature monitor, temperature monitor is situated underneath the pressure pot (perforated metal plate)

or

- Temperature monitor is defective and has to be exchanged

- Control compressed air connection
- Turn on the timer
- Magnetic valve has to be exchanged by an expert

- Clean and grease the lid O-ring
- Renew the lid O-ring

- Reduce the input pressure

- Unit has to be cooled down

## 2. INSTRUCTION

- |  |  |   |
|--|--|---|
| ▪ Unit does not cool down automatically      | ▪ Temperature monitor type: 1822 L 20-4 98 °C (temperature monitor) is defective | ▪ Exchange the temperature monitor                |
| ▪ Unit does not keep the holding temperature | ▪ Magnetic valve is defective  | ▪ Magnetic valve has to be exchanged by an expert |
|  | ▪ Magnetic valve is defective  | ▪ Magnetic valve has to be exchanged by an expert |
|  | ▪ Temperature monitor type: 1822 L 20-4 98 °C (temperature monitor) is defective | ▪ Exchange the temperature monitor                |



### 2.14 Spare part list

| No. | Description                         | Item-no.    |
|-----|-------------------------------------|-------------|
| 1   | Handle                              | D-50275     |
| 2   | Lid                                 | D-50479KPL  |
| 3   | Electronics/digital displays        | 55200       |
| 4   | Counter bearing                     | D-50417KPL  |
| 5   | Seal housing                        | 50617       |
| 6   | Housing silver                      | 52571silber |
| 7   | O-ring 164 x 4 mm                   | 51031       |
| 8   | Swinging screw connection manometer | 50170       |
| 9   | Swinging screw connection pressure  | 50170       |
| 10  | Polymerisation pot PM 1             | D-55552     |
| 11  | Temperature monitor                 | 55743       |
| 13  | Bottom plate PM 1-3                 | 510964      |
| 14  | Main switch                         | 50848       |
| 15  | Mains supply/fuse carrier           | 51208       |
| 16  | Pressure control valve              | 50792       |


## 2. INSTRUCTION

| No.   | Description                                 | Item-no.      |
|-------|---|---------------|
| 17    | Manometer                                   | 50702         |
| 18    | Rotary encoder                              | 55200P        |
| 19    | Lid block                                   | 50305         |
| 20    | Spacer lid bolt                             | 51331         |
| 21    | Rear panel of unit PM 1, 1/120 °C           | 510977        |
| 22    | Magnetic valve                              | 50741-2       |
| 23    | Drain tap                                   | 50744         |
| 24    | Pressure reduction valve                    | 50726         |
| 25    | Compressed air supply                       | 50268         |
| 26    | Ventilation valve                           |               |
|       | Power cord EU                               | 51021         |
|       | Power cord UK                               | 51021/51213   |
|       | Power cord USA/JP                           | 51163/51207   |
|       | Heating cartridge                           | 51054         |
|       | Drain hose                                  | 2004          |
|       | Adhesive display foil PM 1+3                | 510209        |
|       | Rectifier magnetic valve                    | 50741B-Brücke |
| 1a    | Handle                                      | D-50319       |
| 12/22 | Magnetic valve                              | 50741         |
| 7a    | O-ring 164 x 4/120 °C                       | 51168         |
| 11a   | Temperature monitor lid                     | 50701         |
| 11b   | Temperature monitor maintaining temperature | 50935         |
| 27    | Signal lamp green                           | 50843         |


### 3. EXPLOSION SKETCH

## 4. OPERATION DIAGRAM

## 5. PNEUMATIC PLAN

## 6. DECLARATION OF CONFORMITY

Low Voltage Directive (2006/95/ EG).

Herewith we declare that the machine described below in its conception and design and in the shape delivered by us is in accordance with the fundamental requirements for safety and health as prescribed in the applicable EC directives .

In the case of any change or modification of the machine, not authorised by us, this declaration becomes invalid.

### **Denomination of the machine:**

CeraPower 2.0

D-34281I-A

### **Type of machine:**

Pressure polymerisation unit

### **Applicable EC directives :**

2006/95/EG Low Voltage Directive

2004/108/EG EMV Directive

### **Applied harmonised standards :**

#### **2006/95/EG Low Voltage Directive**

**EN 61010-1:2010:** Safety requirements for electrical equipment for measurement, control and laboratory use - Part 1:General requirements (IEC 61010-1 2010 +Cor :2011); German version EN 61010-1:2010

**EN 61010-2-010:2003:** Safety requirements for electrical equipment for measurement, control, and laboratory use - Part 2-010: Particular requirements for laboratory equipment for the heating of materials (IEC 61010-2-010:2003); German version EN 61010-2-010:2003

### **2004/108/EG EMV Directive**

**EN 61000-3-2:2010-03 + Corrigendum to EN 61000-3-2.2011 06:** Electromagnetic compatibility (EMC)-Part 3-2: Limits - Limits for harmonic current emissions (equipment input current <= 16 A per phase) (IEC 61000-3-2:2005 + A1 2008 + A2:2009); German version EN 61000-3-2:2006 +A1:2009 + A2:2009

**EN 61000-3-3:2013:** Electromagnetic compatibility (EMC)-Part 3-3: Limits - Limitation of voltage changes, voltage fluctuations and flicker in public low-voltage supply systems, for equipment with rated current <=16 A per phase and not subject to conditional connection (IEC 61000-3-3:2013); German version EN 61000-3-3:2013

**EN 61326-1:2013:** Electrical equipment for measurement, control and laboratory use - EMC requirements- Part 1: General requirements (IEC 61326-1:2012); German version EN 61326-1:2013 Emission accordmg to living area, business and industrial undertakings as well as small enterprises interference stability according to industrial area.

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**Invicon chemical solutions GmbH**