TwinLux

invicon
chemical solutions

Instructions for use
Informations concernant l'utilisation
Gebrauchsinformation

Date information prepared:
10/2007

Manufacturer

Invicon chemical solutions GmbH

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A 6830 Rankweil

Austria

Phone: +43 (5522) 45 3 01-0

Fax: +43 (5522) 45 3 01-10

office@invicon.at

www.invicon.at

Table of Contents

1	Safety	3	5	Maintenance and cleaning	11
1.1	General information	3			
1.2	Signs and symbols	3	6	Troubleshooting	11
1.3	Safety notes	4			
1.4	Intended use	5	7	Technical data	12
1.5	Standards and classifications	5	7.1	Light unit	12
1.6	Labels on the device	5	7.2	Battery	12
			7.3	Power supply	12
			7.4	Ambient conditions	12
			7.5	Storage conditions	12
			7.6	Weight	12
2	Product description	6	8	Delivery range and spare parts	12
3	First time start-up	8	9	Warranty	13
3.1	TwinLux	8	9.1	Warranty terms and conditions	13
3.2	Charging the battery	9	9.2	What is not covered by the warranty	13
4	Operation	10	9.3	Procedure in case of a claim under warranty	13
4.1	Carrying out light-curing cycles	10			
4.2	Extending the curing time during the curing cycle	11	EC Declaration of Conformity		14
			10	Guarantee certificate information	15

1 Safety

1.1 General information

Please read the user manual for your device carefully before using it for the first time. It contains important information on safety and using your device and its components.

1.2 Labels in the user manual

Special warning information can be found in the user manual before dangerous tasks. This information is broken down into the following danger classes:

Danger!

Dangers which could lead to serious injury or death.

Warning!

Dangers which could lead to serious injury or significant property damage.

Caution!

Dangers which could lead to injury or property damage.

In addition to these warnings, there are also general notes containing useful information.

Note!

Contents of the note.

1.3 Safety information

Danger!

Light from the LED light source

Risk of explosion in areas with highly flammable air conditions.

-> Do not operate the device in areas where there is a risk of explosion.

Danger!

Electromagnetic disturbances

May induce malfunctions in pacemakers or other implanted electronic devices.

-> Do not use the device in the vicinity of a pacemaker or any other implanted electronic device.

Warning!

Light from the LED light source

- Do not look directly into the LED light source.
- Only use the device when the protective shield is attached and functioning correctly.
- Always ensure the protective shield is attached correctly.
- Never point the beam of the LED light source directly at eyes.
- Always aim the beam directly at the polymerisation surface.

Warning!

Defective device

Risk of injury due to parts of the device not functioning correctly.

Check the following points before each use:

- The device and its components are in good working order.
- The protective shield is attached correctly.
- The ventilation slots are not obstructed.

Caution!

Exposure of the device to water

Damage to the device due to water exposure.

- Do not hold the device under running water or immerse it in water.
- Keep the device out of the reach of children.
- Only use genuine accessories supplied by the manufacturer.
- In the event of malfunction, request assistance from an authorised customer support service or from Invicon chemical solutions immediately.

1.4 Intended use

The device is intended for the polymerization of light-curing materials for manufacturing and repair in the craft or industrial sectors.

- The device may only be used by trained specialist personnel.
- Always operate the device with light protection (protective shield).

> We recommend wearing protective goggles.
(Optional available)

1.5 Standards and classifications

The device complies with the following standards and classifications:

EC standard: EN 6060-1

Classification according to IEC 60825-1: 2001: LED Product Class 2

Classification according to CIE S009: Risk Group 2

1.6 Marking of the device

The following device characteristics are located on the device and on the power supply:

- Model
- Serial number
- Supply voltage
- Power consumption

Please always provide us with this information when contacting an authorized service center or directly to Invicon chemical solutions. These data enable our customer service to ensure fast and efficient problem resolution.

2 Product description

NOTE!

Light unit shown on the pictures of the Operating Instructions

The pictures below refer to TwinLux unit.

TwinLux with accessories

- | | |
|---------------------|---------------------------|
| 1 Power pack plug | 5 LED light source |
| 2 Electrical socket | 6 Stand |
| 3 Power pack | 7 On/off and timer switch |
| 4 Cooling vents | 8 Light protection shield |

Power pack (3)

TwinLux is powered by means of the power pack (3) supplied together with the light unit. To connect the power pack to the light unit, the power pack plug (1) is connected to the electrical socket (2) at the rear of the light unit.

TwinLux is powered with a battery, which is charged by means of the power pack (3) supplied together with the light unit.

Air vents (4)

TwinLux has been designed for intensive use and continuous operation. The light unit is equipped with overheating protection including air vents (4) to protect it from overheating. The overheating protection system of TwinLux automatically switches on if the temperature of the heat dissipator rises above 40°C and additional air ventilation cools the light unit. If the temperature goes back down below 40°C, the overheating protection system automatically switches off.

LED light source (5)

WARNING!

Light emission of LED source

Direct exposure to the light of the LED light source may cause damage to the eyes.

- > Do not look directly into the LED light source.
- > Use the light unit only with a properly mounted and fully functioning light protection shield.
- > Always make sure that the light protection shield is mounted properly.
- > Do not point the light beam of the LED source directly into the eyes.
- > Always direct the light beam of the LED light source directly towards the surface to be polymerized.

The LED light source (5) comprises two LEDs, which are positioned behind a protective plexiglass window.

Stand (6)

The integrated stand (6) is used to position the light source at an appropriate distance to the polymerization surface during the curing procedure.

On/off and timer switch (7)

The LED light source (5) is switched on and off by means of the on/off and timer switch (7). The individual curing times are chosen by applying varying degrees of pressure to the on/off and timer switch. As soon as the on/off switch is pressed down and the LED source starts emitting light, the curing time indicator (11) lights up or starts blinking.

Light protection shield (8)

The light protection shield prevents uncontrolled light emission from the LED source (5).

LED indicators

- 9 Power on indicator
- 10 Battery charge indicator
- 11 Curing time indicator

Power on indicator (9)

The red light emitting LED (9) on the rear of the light shows whether the unit is connected to the power supply:

Power supply	TwinLux
LED is off	No power supply
Continuously illuminated	Battery is being charged with power

Battery charge indicator (10)

The red light emitting LED (10) on top of the light shows the amount of power left in the battery (see also Section 3.2).

NOTE!

Exact indication of the amount of power left

You are given a more precise indication of the amount of power left in the battery after having carried out two to three curing cycles.

Curing time indicator (11)

The green light emitting LED (11) on top of the unit shows which curing setting has been selected. The unit switches off automatically after the curing time has elapsed (except for continuous operation).

Curing time indicator (11)	TwinLux
Fast blinking	30 seconds
Slow blinking	60 seconds
Continuously illuminated	Continuous light

3 First time start-up

3.1 TwinLux

Checking the battery charge

> Check the amount of power left in the battery.

NOTE!

Accuracy of the remaining time available

The remaining time indications are only valid at room temperatures (approx. 18 – 25 °C) and after the device has been used in continuous operation for approximately five minutes.

Indication on the light unit	Battery charge	Time remaining on the battery
Battery charge indicator (10) blinks 1 time	Full	Approx. 30 minutes left
Battery charge indicator (10) blinks 2 to 3 times	2/3 full	Approx. 20 minutes left
Battery charge indicator (10) blinks 4 to 5 times	1/3 full	Approx. 15 minutes left
Battery charge indicator (10) blinks continuously	Battery low	Approx. 10 minutes left There is enough power left in the battery to carry out the following light-curing cycles: > 10 to 20 curing cycles of 30 seconds each > 5 to 10 curing cycles of 60 seconds each
The light does not switch on when the on/off switch is pressed	Battery empty	Connect the power pack to the light unit and the mains and charge the battery for 7 to 12 hours
The light switches off automatically	Battery defective	Have the light unit checked by an authorized service centre or directly by Invicon chemical solutions

3.2 Charging the battery

NOTE!

Prolonging the life of the battery

Carrying out a full charge and discharge cycle prolongs the life of the battery.

NOTE!

Battery charge level at delivery

The battery is not fully charged upon delivery.

- > Charge the battery for 12 to 14 hours before using it for the first time.
- > Always re-charge the battery to its full capacity (charging time: 7 to 12 hours) after it has been discharged.

- > Place the light unit on a flat, horizontal surface.
- > Connect the plug (1) of the power pack to electrical socket (2) of the light unit.
- > Connect the power pack (3) to the mains.

4 Operation

WARNING!

Light emission of LED source

Direct exposure to the light of the LED light source may cause damage to the eyes.

- > Do not look directly into the LED light source.
- > Use the light unit only with a properly mounted and fully functioning light protection shield.
- > Always make sure that the light protection shield is mounted properly.
- > Do not point the light beam of the LED source directly into the eyes.
- > Always direct the light beam of the LED light source towards the surface to be polymerized.

WARNING!

Defective device

Malfunctioning devices or components may cause injury.

- > Before each use, check the device for the following:
 - The device and its accompanying components are in proper working condition.
 - The light protection shield is correctly mounted on the light unit.
 - The air vents are free of obstruction.

4.1 Carrying out light-curing cycles

Light emission

- 5 LED light source
- 7 On/off and timer switch
- 8 Light protection shield

TwinLux

- > Mount the light protection shield (8) over the LED light source (5).
- > Select the desired curing time by pressing the on/off and timer switch (7) as follows:

On/off switch (7)	Curing time	Curing time indicator (11)
Short touch (< 1 sec)	30 seconds	Fast blinking
Short pressure (approx. 2 sec)	60 seconds	Slow blinking
Sustained pressure (> 2 sec)	Continuous light emission	Continuously illuminated

- > To interrupt or stop the curing cycle, press again on the on/off switch (7).

4.2 Extending the curing time during the curing cycle

The curing time can be extended to the next longer cycle during operation as follows

TwinLux

- > To increase the curing time from 30 to 60 seconds, press the on/off switch (7) for 2 seconds, while the 30-second cycle is in operation.
- > To increase the illumination time from 60 seconds to continuous operation, press the on/off switch (7) for more than 2 seconds, while the 60-second cycle is in operation.

5 Maintenance and cleaning

If the surface of the curing light or LED light source (5) has become soiled, proceed as follows:

- > Switch off the light unit.
- > Disconnect the plug (1) of the power pack from the light unit.

CAUTION!

Cleaning the surfaces of the device

Exposure to water may damage the light unit

- > Do not hold the light unit under running water. Do not dip the light unit into water.
 - > Do not allow fluid to enter the air vents (4).
 - > Allow the device to dry for at least 5 minutes after cleaning.
- > Clean surfaces with a damp cloth (methylated spirits containing detergent).

If the air vents (4) are heavily contaminated with dirt, proceed as follows:

- > Switch off the light unit.
- > Disconnect the plug (1) of the power pack from the light unit.
- > Clean the air vents (4) with compressed air.

6 Troubleshooting

Malfunction	Possible problem	Solution
Light unit switches off prematurely	Battery is defective	Have the light unit checked by an authorized service centre or directly by Invicon chemical solutions

7 Technical data

7.1 Light unit

Light source:	2 LEDs
Wavelength:	440nm – 480nm
Emission duration:	Unrestricted

7.2 Battery

First time charging:	12 – 14 hours
Re-charging:	7 – 12 hours
Cells:	TwinLux: 6, m3: 4
Type:	NiMH, HHR-210AAB2B
Voltage:	TwinLux: 7.2V, m3: 4.8V
Power pack:	12.6A

7.3 Power supply

Power pack:	90 – 260 VAC / 50 – 60 Hz
-------------	---------------------------

7.4 Ambient conditions

Temperature:	+10°C – +40°C
Relative air humidity:	30% – 85%

7.5 Storage conditions

Temperature:	-10°C to +70°C
Relative air humidity:	10% – 90%
Air pressure:	500hPa – 1060hPa

7.6 Weight

TwinLux	394 g (without power pack)
---------	----------------------------

8 Delivery range and spare parts

- 1 x TwinLux Unit
- 1 x Power Pack
- 1 x Light Protection Shield

- 1 x Operating Instructions
- 1 x Warranty Form

TwinLux

Power Pack

Light Protection Shield

9 Warranty

9.1 Warranty terms and conditions

Invicon chemical solutions gives a warranty of 2 (two) years on all TwinLux units as of the date of purchase. This also applies to light units purchased from an Invicon retailer or importer. The warranty is limited to defects related to the manufacture and material of the product.

Invicon, in its sole discretion, may repair or replace the defective parts of the product free of charge within the warranty period.

In general, the „General Terms and Conditions of Sale and Delivery“ shall apply.

This warranty does not provide for the replacement of entire TwinLux.

9.2 What is not covered by the warranty

The following will void this warranty:

- The apparatus and its components were not used in compliance with the safety instructions and directions provided in these Operating Instructions.
- The apparatus and its components were not used for the intended purpose.
- The apparatus or its components were altered by the user or other third parties.
- Assembly, settings, alterations, additions and repairs were not carried out by specialized Invicon personnel or authorized specialist.
- The electrical facilities of the room in which the apparatus was used did not comply with the current local laws and regulations.
- The conditions in the premises where the apparatus was used, kept or stored did not comply with the conditions listed in the Technical Data section of these Operating Instructions.
- Damage arising from catastrophic events, such as impact, accident, vandalism and events beyond the control of the manufacturer.

The following items are not covered by this warranty:

- LED light source
- LED indicators
- Switch
- Accompanying components/accessories

The following damage is not covered by the warranty:

- Damage caused during transportation
- Damage arising from improper use or negligence
- Damage arising from wrong usage of electrical supply and voltage

9.3 Procedure in case of a claim under warranty

The warranty is only valid if the warranty form supplied with the apparatus is fully filled in and returned together with the receipt (including the date of purchase) to Invicon chemical solutions or to an authorized Invicon dealer within fifteen (15) days as of the date of purchase.

In case of a claim under warranty, the customer shall put the light unit at the disposal of Invicon or the company's distribution partner and provide the following information:

- Name, address, telephone number, email of the owner
- Name and address of the Invicon distribution partner from whom the light unit was purchased
- Serial number of the light unit
- Receipt (delivery note)
- Detailed description of the malfunction

The data and descriptions contained in these Operating Instructions are not binding and may be varied by the manufacturer at any time without previous notice.

EC Declaration of Conformity

Invicon chemical solutions GmbH
Schweizer Strasse 96
6830 Rankweil
Austria


We, the undersigned, declare under our sole responsibility that the devices specified below

TwinLux

Has been developed, built and manufactured in compliance with the following EC Directives (as well as in compliance with the provisions of domestic law adopted in the field governed by these Directives):

- EN 61010-1
- EN 50081-2
- EN 50082-2

- 73/23 EWG
- 89/336/EWG

Rankweil, 1 October 2007

Gerhard Zanghellini, General Manager

Guarantee certificate information

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