TwìnLux Dual

Instructions for use

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office@invicon.at
www.invicon.at

Contents

1 Safety	3	7 Technical data	15
1.1 General instructions	3	7.1 Light sources	15
1.2 Labels in the user manual	3	7.2 Battery	15
1.3 Safety information	4	7.3 Power supply	15
1.4 Intended use	5	7.4 Weight	15
1.5 Standards and classifications	5	7.5 Environment / Storage conditions	15
1.6 Labels on the device	5	7.6 Risk Group	15
2 Device description	6	8 Scope of delivery and spare parts	16
3 Start-up operation	9	9 Guarantee	16
3.1 Checking the charge level	9	9.1 Guarantee conditions	16
3.2 Charging the battery	11	9.2 Guarantee exclusion	16
4 Operation	11	9.3 Guarantee claim procedure	17
4.1 Applying light	12	EC Declaration of Conformity	18
4.2 Increasing the exposure time during exposure	13	10 Guarantee certificate information	19
5 Maintenance and cleaning	14		
6 Troubleshooting	14		

1| Safety

1.1 General information

Please read the user manual for your device carefully before using it for the first time. It contains important information on safety and using your device and its components.

1.2 Labels in the user manual

Special warning information can be found in the user manual before dangerous tasks. This information is broken down into the following danger classes:

Danger!

Dangers which could lead to serious injury or death.

Warning!

Dangers which could lead to serious injury or significant property damage.

Caution!

Dangers which could lead to injury or property damage.

In addition to these warnings, there are also general notes containing useful information.

Note!

Contents of the note.

1.3 Safety information

Danger!

Light from the LED light source

Risk of explosion in areas with highly flammable air conditions.
-> Do not operate the device in areas where there is a risk of explosion.

Danger!

Electromagnetic disturbances

May induce malfunctions in pacemakers or other implanted electronic devices.
-> Do not use the device in the vicinity of a pacemaker or any other implanted electronic device.

Warning!

Light from the LED light source

- Do not look directly into the LED light source.
- Only use the device when the protective shield is attached and functioning correctly.
- Always ensure the protective shield is attached correctly.
- Never point the beam of the LED light source directly at eyes.
- Always aim the beam directly at the polymerisation surface.

Warning!

Defective device

Risk of injury due to parts of the device not functioning correctly.

Check the following points before each use:

- The device and its components are in good working order.
- The protective shield is attached correctly.
- The ventilation slots are not obstructed.

Caution!

Exposure of the device to water

Damage to the device due to water exposure.

- Do not hold the device under running water or immerse it in water.
- Keep the device out of the reach of children.
- Only use genuine accessories supplied by the manufacturer.
- In the event of malfunction, request assistance from an authorised customer support service or from Invicon chemical solutions immediately.

1.4 Intended use

The device is intended to polymerise light-curing materials of HyCeram® Bond II or Bond THX NF, LuxArt® and Luxon® design materials.

- The device may only be used by trained specialists.
- Always use light protection (the protective shield) when operating the device.

> We recommend wearing protective goggles.
(Optional available)

1.5 Standards and classifications

EC standard: EN62471

IEC classification Risk group 2

CIE classification Risk group 2

1.6 Labels on the device

The following device information can be found on the device itself:

- Model
- Serial number
- Supply voltage
- Power consumption

Please always provide this information if you contact an authorised customer support service or Invicon chemical solutions directly. This data allows us to ensure that problems are solved quickly and efficiently.

2| Device description

The TwinLux Dual is a mains-independent, portable device with a built-in battery. It is used to cure the following light-curing materials:

- HyCeram® Bond II or Bond THX NF materials
- LuxArt® design composite
- Luxon® design acrylate

More detailed descriptions of HyCeram®, LuxArt® and Luxon® can be found in the corresponding instructions for use.

Two high-efficiency Royal Blue LED diodes with a wavelength range of 425 - 470 nm (peak 450 nm) are used as light sources.

One UV LED diode with a wavelength range of 350 - 375 nm (peak 365 nm) is also used.

The total radiant power is 1.5 W.

This means that the entire LED light spectrum emitted by TwinLux Dual is used to full effect in the activation of the HyCeram® Bond II or Bond THX NF, LuxArt® or Luxon® initiators. This produces excellent polymerisation results with low energy consumption and minimal heat loss.

The LED light sources can be switched on and off using the multi-function button (4).

Three different preset exposure times can be selected, depending on how the button is pressed.

The following exposure times can be chosen:

- 30 seconds
- 60 seconds
- Continuous exposure

As soon as the function button is pressed and the light is active, the exposure indicator will light up or start to blink.

The exposure times can be activated as follows:

Multifunction button action (7)	Exposure time	Exposure time indicator (11)
Short tap	30 seconds	Fast blinking
Short press	60 seconds	Slow blinking
Long press	Continuous exposure	Steady light

Note!

Illustration in the user manual

The following illustration always refers to the TwinLux Dual device.

TwinLux Dual with accessories:

1 Power supply unit plug
2 Power connector on device

3 Mains connection

4 Overheat protection with ventilation slots

5 LED light source

6 Stand

7 Function button on/off

8 Protective shield

Mains connection (3)

The mains connection supplies power to the TwinLux Dual and charges the battery.

Overheat protection (4)

The TwinLux Dual is designed for intensive, continuous use. The device features overheat protection with ventilation slots (4) in order to prevent overheating.

LED light source (5)

Warning!

Light from the LED light source.

Risk of eye injury from direct exposure to light.

- Do not look directly into the LED light source.
- Only use the device when the protective shield is attached and functioning correctly.
- Always ensure the protective shield is attached correctly.
- Never point the beam of the LED light source directly at eyes.
- Always aim the beam directly at the polymerisation surface.

The LED light source (5) contains two Royal Blue LEDs and one UV LED behind a protective Plexiglas pane.

Royal Blue LED: Wavelength range of 425 - 470 nm (peak 450 nm).

UV LED: Wavelength range of 350 - 375 nm (peak 365 nm).

Stand (6)

The stand (6) is used to secure the unit.

Function button (7)

The multifunction button (7) is used to switch the LED light source (5) on and off. As soon as the function button is pressed and the light is active, the exposure time indicator (11) will light up or start to blink.

Protective shield (8)

Attaching the protective shield prevents uncontrolled light emission from the LED light source (5).

LED indicators

- 9 Power supply indicator
- 10 Charge level indicator
- 11 Exposure time indicator

Power supply indicator (9)

The red LED (9) on the back of the device shows the power supply status:

Power supply indicator

LED off
Steady light

TwinLux Dual

No power supply
Battery is charging

Charge level indicator (10)

The red LED (10) on the top of the device indicates the battery charge level of the light unit (see also 3.2).

Note!

Accurate charge level indication

The charge level will be indicated more accurately if you carry out 2 - 3 exposures beforehand.

Exposure time indicator (11)

The green LED (11) on the top of the device indicates the selected exposure time. The device will automatically switch off after the selected exposure time is up (unless continuous exposure is selected).

Multifunction button action (7)	Exposure time	Exposure time indicator (11)
Short tap	30 seconds	Fast blinking
Short press	60 seconds	Slow blinking
Long press	Continuous exposure	Steady light

3| Start-up operation

- > Place the device on a horizontal surface.
- > Connect the power supply unit plug (1) to the power connector.

3.1 Checking the charge level

- > Check the battery charge level.

Note!

Validity of time indications

The time indications are only accurate at room temperature (around 18 – 25° C) and after around five minutes of continuous use.

Device status	Battery status	Battery capacity
Charge level indicator (10) blinks 1 x	Full	Around 30 minutes of use remaining
Charge level indicator (10) blinks 2 - 3 x	2/3 full	Around 20 minutes of use remaining
Charge level indicator (10) blinks 4 - 5 x	1/3 full	Around 15 minutes of use remaining
Charge level indicator (10) blinks continuously	Battery nearly empty	Around 10 minutes of use remaining The following exposure cycles are still possible: >10 - 20 exposure cycles of 30 seconds each >5 - 10 exposure cycles of 60 seconds each
Light does not activate when button is pressed	Battery drained	Connect the power supply unit and charge the battery for 7 to 12 hours
Device turns off automatically	Battery defective	Have the device examined by an authorised specialist or directly by Invicon chemical solutions.

3.2 Charging the battery

Note!

Increased battery lifespan

Charging the battery fully and then discharging it fully will increase the lifespan of the battery.

Note!

Battery charge level on delivery

The battery is not fully charged on delivery.

- Fully charge the battery for 12 - 14 hours before first use.
- Always fully charge the battery for 7 - 12 hours after the first discharge.

4| Operation

Warning!

Light from the LED light source

Risk of eye injury from direct exposure to light.

- Do not look directly into the LED light source.
- Only use the device when the protective shield is attached and functioning correctly.
- Always ensure the protective shield is attached correctly.
- Never point the beam of the LED light source directly at eyes.
- Always aim the beam directly at the polymerisation surface.

Warning!

Defective device

Risk of injury due to parts of the device not functioning correctly.

Check the following points before each use:

- The device and its components are in good working order.
- The protective shield is attached correctly.
- The ventilation slots are not obstructed.

4.1 Applying light

Exposure

- 5 LED light source
- 7 Function button
- 8 Protective shield

TwinLux dual

Attach the protective shield (8) over the LED light source (5).
Press the function button (7).

Multifunction button action (7)

Short tap
Short press
Long press

Exposure time

30 seconds
60 seconds
Continuous exposure

Exposure time indicator (11)

Fast blinking
Slow blinking
Steady light

Once light exposure has started, briefly press the function button (7) again to stop the exposure.

Please observe the following exposure times and distances

Material	Exposure time	Distance
HyCeram® Bond II or Bond THX NF	120 seconds	2 cm
LuxArt® Design	120 seconds	2 cm
Luxon® Design	120 seconds	2 cm

Note!

Ensure that all material surfaces are exposed to light carefully.
(Exposure Time, distance).

4.2 Increasing the exposure time during exposure

It is also possible to increase the exposure time by one level during exposure:

- To increase the exposure time to 60 seconds, when an exposure time of 30 seconds has been chosen, hold down the function button (7) for two seconds.
- To increase the exposure time to continuous exposure, when an exposure time of 60 seconds has been chosen, hold down the function button (7) for longer than two seconds.

5| Maintenance and cleaning

If the surface of the device or the LED light source (5) gets dirty:

- > Switch off device
- > Disconnect the power supply unit plug (1)

Caution!

Cleaning the surface of the device

Damage to the device due to water exposure.

- Do not hold the device under running water or immerse it in water.
- Do not let liquids enter the ventilation slots (4).
- Leave the device to dry for at least five minutes after cleaning.
- Clean the surface with a damp cloth (alcohol-based cleaner).

If the ventilation slots (4) are very dirty:

- Switch off the device.
- Disconnect the power supply unit plug (1).
- Clean ventilation slots (4) with compressed air.

6| Troubleshooting

Problem	Potential cause	Solution
Device turns off automatically	Battery	Have the device examined by an authorised specialist or directly by Invicon chemical solutions.
Device turns off automatically	Malfunction	Have the device examined by an authorised specialist or directly by Invicon chemical solutions.

7| Technical data

7.1 Light sources

LEDs	2 x LED Royal Blue 425 - 470 nm (peak 450 nm) 1 x LED UV 350 - 375 nm (peak 365 nm)
Total radiant power	1,5 W

7.2 Battery

Initial charge time	12 - 14 hours
Recharge time	7 - 12 hours
Number of cells	4
Type	NiMH, HHR, 210AAB2B
Voltage	4.8 V
Power supply unit	DC 12V/1A

7.3 Power supply

Power supply unit	90 - 260 VAC / 50 - 60 Hz
-------------------	---------------------------

7.4 **Weight:** 394 g (not including power supply unit)

7.5 Environment/Storage conditions

Temperature	-10° C - +70° C
Rel. humidity	10 % - 90 %
Air pressure	500 hPa - 1060 hPa

7.6 **Risk group:** 2 (Classification in accordance with DIN EN 62471 300 nm - 1400 nm)

8| Scope of delivery and spare parts

1x TwinLux Dual light unit
1x power supply unit
1x protective shield
1x user manual
1x guarantee certificate

9| Guarantee

9.1 Guarantee conditions

Invicon chemical solutions guarantees all TwinLux Dual light units sold (including by an Invicon reseller or importer) against material or production faults, effective for two years from the date of purchase. During the guarantee period, Invicon chemical solutions agrees to repair parts it considers defective free of charge or to replace them at its own discretion. In general, the relevant valid version of the General Terms and Conditions of Business and Delivery shall apply.

Replacement of the entire TwinLux Dual light unit is not possible.

9.2 Guarantee exclusion

In the following cases, guarantee claims shall be void:

- The device or the associated components were not used in accordance with the safety rules and instructions described in these instructions for use.
- The device or the associated components were not used in accordance with the intended use.
- The device or components were structurally modified without authorisation.

- Assembly, settings, modifications, extensions or repairs were not carried out by Invicon specialist personnel or authorised specialist companies.
- The electrical system in the rooms in which the device is operated does not correspond to the valid laws or rules.
- The ambient conditions in which the device is stored or used do not correspond to the conditions listed in the “Technical data” section.
- Disaster situations as a result of foul play, accident, vandalism or force majeure.

The following components are excluded from the guarantee:

- LED light sources
- LED indicators
- Switches
- Accessories

The following damage is excluded from the guarantee:

- Transport damage
- Damage resulting from improper use or neglect
- Damage resulting from connection to a voltage that the device is not designed for

9.3 Guarantee claim procedure

The guarantee is only valid if the guarantee certificate provided with the unit is completed in full, including proof of purchase with date of purchase, and delivered to our company headquarters or an authorised Invicon chemical solutions sales partner within 15 days of the date of purchase.

In the event of a guarantee claim, the customer must provide the device to Invicon chemical solutions or the Invicon sales partner, and must supply the following information:

- Name, address, phone number and e-mail address of the buyer
- Name and contact details of the Invicon distributor
- Serial number of the device
- Proof of purchase for the device (delivery note)
- Detailed description of the malfunction

The information contained in this description is non-binding and may be changed at any time without prior notice.

EC Declaration of Conformity

pursuant to EU Low Voltage Directive 2014/35/EU, dated 26/2/2014

Manufacturer: oehri electronic AG | Industriering 35 | 9491 Ruggell, Liechtenstein

Description and identification of interchangeable equipment:

Product: TwinLux-Dual
Task: NA16-002406 and following
Function: Blue/UV portable exposure unit

It is expressly declared that the interchangeable equipment complies with all relevant conditions of the following EU directives:

2014/35/EU: (Low Voltage Directive) Directive of the European Parliament and of the Council of 26 February 2014 on the harmonisation of the laws of the Member States relating to the making available on the market of electrical equipment designed for use within certain voltage limits

2011/65/EC RoHS: Restriction of the use of certain hazardous substances

Reference of applied harmonised standards:

EN61000-6-1:2007 Electromagnetic compatibility (EMC) - Part 6.1 Generic standards
EN61000-6-3:2006-3 Emission standard and immunity for residential, commercial and industrial environments
EN 62471 300 to 1400 nm, risk group 2 (moderate risk):
Labelling is required.

EC Declaration of Conformity

pursuant to EU Low Voltage Directive 2014/35/EU, dated 26/2/2014

oehri electronic AG
Ruggell, 22 November 2016

Manfred Laim
Managing Director

Guarantee certificate information

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Invicon chemical solutions GmbH